



Ajax Dental Supplies Pty Ltd.

Headquarter ADD: 10/19 Cornhill Street Ferntree Gully VIC 3156 (Showroom / Service Centre)
NSW Branch: Unit 18 25-33 Alfred Road, Chipping Norton NSW 2170 (Showroom / Service Centre)
Tel: 03 9879 3801 Mobile: 04 0233 9905 After sales hotline: 0478 635 358

E-mail: info@ajaxdent.com.au



<http://www.ajaxdent.com.au>



Ajax Dental


AJ18 Dental Unit The Pursuit Of Excellence

www.ajaxdent.com.au


E-Mail: info@ajaxdent.com.au




Fast Motor
The backrest can be fully restored in 9 seconds.




Cast aluminum alloy patient chair
The integrally cast metal chair frame makes the system sturdy and durable.




**Dual water purification system
Dual air purification system**




Rotatable spittoon
The spittoon is rotatable, convenient for patient to rinse.




Doctor instrument
·Luxury suspension system with anti-backflow system (15 functions)
·Chair-instrument mutual lock system.
·Built-in X-ray viewer
·3 Air tubing hose
·1 Three way syringe



Dental light
Whale LED light with touch-free sensor and two operating models.



Multi-function foot control
Including chair position control, instrument panel control, Rinse and cup filling control.



Water/Air/Power one switch on/off control system
After finishing daily work, the doctor should turn off the switch to prevent excessive water pressure at night and prolong the life of the pipeline.

Standard Configurations			
Dental Light	Whale LED light with touch-free sensor and two operating models	Side Box	Cast aluminum water unit frame with swivel side box
Chair	Chair safety switch		24V water heater
	Cast aluminium backrest stand and base, cast iron chair frame		Rotatable ceramic cuspidor
	3 programmable positions		Imported Solenoid valve system
	Press type articulating dual axis headrest		1500ml Water purification system
	Fiber leather upholstery		1300ml Disinfectant system
	Swing-out armrest		Imported hydrolytic resistance pipeline
	Electric motor system (High speed backrest motor system)		Luxury hanging system with anti backflow system (15 functions)
Assistant Instrument	Suction delay off system	Doctor Instrument	Chair-instrument mutual lock system
	10 function control panel (assistant side)		Built-in X ray viewer
	Double joint assistant stents		Air turbine hose 3pcs (import from US, stainless steel connector)
	1 HVE and 1 SE with adjustable switch		1 Three way syringe
	1 Three way syringe	Junction Box	Dual water purification system
	With assistant tray		Dual air purification system
Foot Pedal	11 functions foot control		Water/Air/Power master switch control system
	Dental Stool		D3 dental stool with 6 functions
			Power box with main control box

