

Configuration

Operating System	Ajax OS M2 (Ajax System)
	Safety switch
	Soft start and stop
Infection Control	One click water line disinfection system
Chair	Fan-shaped backrest
	Fiber Leather upholstery
	Swing-back armrest
Delivery System	Silicon Pad
	Built-In X-ray viewer
	Master air/Power/Water switch
Assistant Unit	Single joint
	Rotatable
Side Box	Side box 1300ml water bottle
	Detachable glass cuspidor
	Hydrolytic resistance water lines
Dental Light	Pebble II
Foot Control	FT-6A Multifunctional Foot Control
Junction Box	Outside junction box with air/water filters and regulators
	Sealed drainage device
Dental Stool	D3 stool